

Babi.Plus® nTest Lung

Neonatal Lung Simulator

The Babi.Plus nTest Lung is an easy to use neonatal test lung for demonstration and testing applications with infant and neonatal ventilator support devices. The nTest Lung consists of two individual silicone lungs and a Y connector. The nTest Lung is the ideal instrument for training and demonstration in addition to performing ventilator verification testing.

Specifications

Compliance	Approximate 2.0 mL/cm H ₂ O
Tidal Volume	25 ml each silicone bag, 50 ml total
Resistance	300 cm H ₂ O/L/s at 1 L/sec
Interface Connector	Y connector with 15/10 mm port
Straight Connector	10 mm OD
Pressure Monitoring Port	4.3 mm ID
Weight	43 grams

Ordering Information

2010-300	Babi.Plus nTest Lung, R300 (1 each)
----------	-------------------------------------

Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit www.respiralogics.com or email us at information@respiralogics.com.



FEATURES

Two lung feature with Y connector to mimic anatomical design

Standard ISO connections

Lightweight, easy to use and portable

Manufactured of high performance material

Packaged in convenient zippered canvas bag for storage and transport



310 Twin Oaks Valley Road • Suite 107-336
San Marcos, CA 92078 • 760 560 3555