

# Jimmy Proof Straps™

## Ventilator circuit security strap to prevent accidental disconnect

Jimmy Proof Straps are an anti-disconnect device intended to secure the ventilator circuit connector to the tracheostomy tube holder to help reduce risk of accidental disconnects. The Straps are made of cotton with reinforced stitching to prevent stretch and rolling of the strap. The ventilator connector threads through the center opening of the Strap and the tapered ends of the Strap attach easily to the tracheostomy tube collar. Engineered to hold its shape with hook tab connectors to ensure quick release when disconnect is indicated.

The Jimmy Proof Straps are intended for use in critical care units, special procedure units as well as alternate care facilities and home care.

Jimmy Proof Straps were developed with special care by a dad wanting to "make a difference" for his son and change the quality of life for all tracheostomy patients. His engineering expertise and commitment to finding a better solution is found in the Jimmy Proof Straps and Danny Ties™.


### Ordering Information

DWL105005	Jimmy Proof Straps (250/case)
DWL105005B	Jimmy Proof Straps (25/box)

### Companion Products

DWL105001	Danny Ties, Small
DWL105002	Danny Ties, Medium
DWL105003	Danny Ties, Large
SP01	Sil.Flex Stoma Pad, Size 1
SP02	Sil.Flex Stoma Pad, Size 2
SP03	Sil.Flex Stoma Pad, Size 3
DP00	Sil.Flex TC Pad, 5 mm
DP01	Sil.Flex TC Pad, 10 mm

Shown with universal adaptor  
(universal adaptor not sold with  
Jimmy Proof Strap)



Respiralogics products are available from Respiralogics and our specialty respiratory care distributors. For our latest product information visit [www.respiralogics.com](http://www.respiralogics.com) or email us at [information@respiralogics.com](mailto:information@respiralogics.com).



**respiralogics™**

310 Twin Oaks Valley Road • Suite 107-336  
San Marcos, CA 92078 • 760 560 3555



### FEATURES

One piece design, constructed of soft cotton with reinforced stitched edges preventing stretching and rolling

Strengthened center hole of the Strap provides a stable opening for ventilator connector

Unique tapered, hook tabs easily adhere to Danny Ties and many commercially available tracheostomy collars

Cotton material absorbs moisture and wicks away from patient

Single patient use item, individually packed in environmentally friendly packaging